

Acute Myocardial Infarction (AMI) Screening 1 Item



Name: Cardiac Troponin I/Creatinekinase-MB/Myoglobin Detection Kit (Lateral Flow Method)

Model: YC-HM2220112 Certificate: ISO13485:2016, EC Declaration of Conformity

Sample Type: Whole blood, Plasma, Serum

Packing Specification: 1 Test/Kit, 20 Tests/Kit

Application: A blood sample is taken without an instrument and quickly assessed for chlamydia pneumoniae antibodies within 15 minutes to determine if you have contracted pneumonia. By taking a blood sample and rapidly detecting and interpreting chlamydia pneumoniae antibodies in a person within 15 minutes.