

Infectious Disease (Qualitative) Virus Screening 12 Items



Name: Hepatitis A Virus Antibody Detection Kit (Lateral Flow Method)
Model: YC-HS2220101 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: A blood sample can be used to quickly screen for HAV antibodies within 15 minutes without an instrument to determine if you have been infected with HAV disease.



Name: Hepatitis B Virus (HBsAg) Surface Antigen Detection Kit (Lateral Flow Method)
Model: YC-HS2220102 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: A blood sample is taken to quickly screen for HBV (HBsAg) within 15 minutes without equipment.



Name: Hepatitis C Virus Antibody Detection Kit (Lateral Flow Method)
Model: YC-HS2220103 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: A blood sample can be taken to quickly screen for hepatitis C antibodies within 15 minutes without an instrument to determine if you have been infected with hepatitis C virus disease.



Name: Hepatitis E Virus Antibody Detection Kit (Lateral Flow Method)
Model: YC-HS2220104 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: Rapid screening for hepatitis E antibodies within 15 minutes through blood sampling, without equipment, to determine whether people have been infected with hepatitis E virus disease.



Name: HAV/HEV Antibody Detection Kit (Lateral Flow Method)
Model: YC-HS2220105 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: Rapid screening for HAV and HEV antibodies within 15 minutes through blood sampling without equipment to determine whether you have been infected with HAV or HEV disease.



Name: Toxoplasma Antibody Detection Kit (Lateral Flow Method)
Model: YC-2220106
Certificate: ISO13485:2016, CE
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit

Application: A blood sample can quickly screen for toxoplasma in 15 minutes without equipment.



Name: Treponema Pallidum Antibody Detection Kit (Lateral Flow Method)
Model: YC-HS2220107 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: Rapid screening for antibodies to Treponema Pallidum within 15 minutes through blood sampling, without equipment, to determine whether you have been infected with Treponema Pallidum disease.



Name: Gonorrhea Antigen Detection Kit (Lateral Flow Method)
Model: YC-HS1550108 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Male urethral swab sample, Female cervical swab sample
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: Rapid screening for Gonorrhea within 15 minutes by means of a urethral swab or cervical swab sampling without equipment can determine infection with Gonorrhea disease.



Name: Monkeypox Virus Antigen Detection Kit (Lateral Flow Method)
Model: YC-HS1660109 **Certificate:** ISO13485:2016, EC Declaration of Conformity
Sample Type: Skin lesion tissue
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: The detection of suspected monkeypox virus infection is determined by rapid screening of monkeypox virus antigens within 15 minutes without equipment through skin lesion tissue sampling.



Name: Monkeypox Virus Antibody Detection Kit (Lateral Flow Method)
Model: YC-HS2220110
Certificate: ISO13485:2016, EC Declaration of Conformity
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: A blood sample is taken to quickly screen for monkeypox antibodies within 15 minutes without equipment to determine if you have been infected with monkeypox virus disease.



Name: HIV1/ 2 Antibody Detection Kit (Lateral Flow Method)
Model: YC-HM2220111 **Certificate:** ISO13485:2016
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: Through blood sampling, without equipment, rapid screening of HIV I, II and O, antibodies within 15 minutes, to determine whether people have been infected with HIV 1 and 2 diseases.



Name: HIV Ag/Ab Detection Kit (Lateral Flow Method)
Model: YC-HM3220212 **Certificate:** ISO13485:2016
Sample Type: Whole blood, Plasma, Serum
Packing Specification: 1 Test/Kit, 20 Tests/Kit
Application: Through blood sampling, without equipment, rapid screening of HIV antigens, antibodies within 15 minutes, to determine whether the infection of HIV disease.